

TabTimer MK3 Careousel Automatic Pill Dispenser



TabTimer Careousel Automatic Pill Dispenser The Pivotell® Automatic Pill Dispenser reminds the user by means of an alarm and flashing light when to take their medicines.



- **1-28 Daily Alarm Times** - Accomodates the most complicated of medication regimes
- **29 Compartments** - For easy pharmacy filling (one compartment is left empty)
- **Dual / Multi Dose programming** - Two or more compartments can be programmed to alarm at the same time. For those with many or large pills
- **Start delay function** - Delays activation to allow pharmacists or family members to easily fill in advance away from the user's home
- **Pharmacy label carrier and label case** - supplied to comply with pharmacy labelling requirements
- **Upside down feature** - The medication tray rotation is disabled if the dispenser is left upside down, to avoid pill spillage. The alarm will still sound at a dosage time.
- **Day / time discs supplied** for 1 / 2 / 3 / 4 times daily dosage, and one blank disc for five or more doses per day
- **Battery life** - up to 6-12 months
- **Battery level meter** and warning indicator
- **Secure battery case cover** to deter tampering
- **Four different alarm sounds**
- **Flashing light** for those with hearing impairment
- **Lockable** - only the medication due can be accessed by the user
- **12/24 Hour Time** - The clock can be set to either 12 or 24 hour display.
- **Automatic daylight savings** - The clock is self adjusting for daylight savings, this can be turned on/off as required

- **Security** - An alarm signal (different to the pre set alarm sound) is activated if any fault should arise which effects the rotation of the dispenser
 - **Testing** - The Pivotell dispenser has an inbuilt test function which should be activated periodically
 - **Lock** - The Pivotell dispenser is equipped with a lock to ensure that only the medication which is due is accessible to the user **Weight** - with batteries 480 gms
 - **The Dispenser Box Includes:**
 - 1x Mark 3/11 PivoTell® Automatic Pill Dispenser
 - 2x time discs for 1-4 times a day regime
 - 2x blank time discs for more than 4 times a day
 - Key and battery cover securing screw
 - Prescription booklet and holder to apply pharmacy labels or medication information
 - 4 x AA batteries
 - **Optional Accessories:**
- Filling aids** to assist with accurate filling of the trays.
- **Spare trays, lids and time discs** for off-site filling by family members or pharmacists
 - **Disposable tray system** for off-site filling by pharmacists
 - **A Pivotell Dispenser Tipper** for those with limited strength or wrist dexterity and those who might forget to turn the dispenser upright after dispensing their medication

| Code | Display Name |
|-----------|---|
| DLG275000 | TabTimer MK3 Careousel Automatic Pill Dispenser - Solid White Lid |

