

Large Capacity Medinizer



e-pill | Large Capacity Weekly Organiser and Medication Reminder A practical and durable weekly medication management system with seven pillboxes and an integrated timer that allows up to six daily alarms. Designed for convenience at home, work, or while traveling — helping you stay on schedule and never miss a dose.



- **Weekly Organiser:** Includes seven individual pillboxes labeled Monday to Sunday for a full week of medication management.
- **Four Daily Compartments:** Each pillbox has four large sections labeled Morning, Noon, Evening, and Night to organise daily doses.
- **Large Capacity:** Each compartment holds up to 25 aspirin-sized tablets, ideal for multiple medications or supplements.
- **Braille Markings:** Each pillbox includes tactile labeling for accessibility and ease of use.
- **Durable Construction:** High-quality design with secure click-on lids and smooth sliding mechanism for easy access.
- **Built-in Timer:** Comes with a medication timer allowing up to six daily alarms to remind you when it's time to take your medication.
- **Easy Alarm Setup:** Simply slide the switches to set alarms on the hour; alarms automatically repeat daily.
- **Audible and Visual Alerts:** At alarm time, a 15-second beeping alarm sounds; if not stopped, it repeats after 5 minutes. A flashing indicator provides a visual alert for 5 minutes.
- **Simple to Use:** Press the stop button after taking your medication to silence the alarm.
- **Secure and Convenient:** Unique click-on lid design keeps pills safely stored while ensuring easy access when needed.
- **Portable Design:** Encased in a stylish vegan leather carry case for discretion

and protection.

- **Organize an entire week's medication with one setup:** The e-pill Large Capacity Weekly Organizer simplifies complex medication schedules, ensures timely reminders, and offers portability with a discreet design.

Specifications:

- **Canvas Bag Size:** 6.5 in × 6.5 in × 2 in
- **Daily Pillbox Size:** 3.5 in × 2 in × 0.75 in
- **Total Weight:** 15 oz
- **Pill Capacity:** 25 aspirin-sized tablets per compartment
- **Color:** Blue canvas bag with clear pillboxes
- **Battery:** 2 × AAA batteries (included)
- **Battery Life:** 6–12 months, depending on the number of daily alarms set

Code	Display Name
DLG275040	Large Capacity Medinizer Weekly Organizer Reminder