

Aspire ComfiGel Cushion



The Aspire ComfiGel Cushion is designed to deliver adaptive comfort, ergonomic support, and all-day freshness. Its contoured shape provides better support, while advanced hyper-elastic gel technology gently conforms to the body to enhance sitting comfort across a wide range of environments at home, in the office, or on the go.



- **Ergonomic Support Design:** Contoured to support the hips and lower back, helping encourage healthy posture, improve sitting balance, and reduce strain during extended periods of sitting.
- **Adaptive Gel Comfort:** Hyper-elastic gel conforms to your natural body shape to provide a personalised sitting experience. The gel structure distributes weight evenly to help support sensitive areas.
- **Easy-care Washable Cover:** Soft, breathable outer cover is fully removable and machine washable for easy cleaning and long-lasting hygiene. Durable stitching and a quality zip provide reliable everyday use.
- **Cool and Breathable Materials:** Open-cell gel combined with a breathable fabric cover promotes airflow, helping to prevent heat and moisture build-up for all-day temperature control and freshness.
- **Lightweight and Portable:** Compact, lightweight construction makes it easy to transfer the cushion between chairs, vehicles, mobility devices, or workspaces, ideal for comfort wherever life takes you.
- **Hyper-elastic Gel:** Two-inch thick gel cushion adapts seamlessly to your body, delivering consistent comfort and support throughout the day.
- **Washable Cover:** Removable, machine-washable cover for simple maintenance and hygiene.
- **Temperature Control:** Breathable materials help prevent heat build-up, keeping you cool during prolonged sitting.
- **Lightweight Design:** Easy to reposition, carry, and use across different seating environments.

Warranty Information

- 1 Year

Code	Display Name	Total Length	Total Width	Total Height
PTC530390	Aspire ComfiGel Cushion - 20" x 18"	20inch	18inch	50mm
PTC530100	Aspire ComfiGel Cushion - 16" x 16"	16inch	16inch	50mm
PTC530200	Aspire ComfiGel Cushion - 16" x 18"	16inch	18inch	50mm
PTC530290	Aspire ComfiGel Cushion - 18" x 16"	18inch	16inch	50mm
PTC530300	Aspire ComfiGel Cushion - 18" x 18"	18inch	18inch	50mm