

Aspire Deluxe Bedside Commode



- Allows for toileting in the bedroom for users unable to travel between the bedroom and bathroom
- Available in Champagne or Slate Vinyl
- Easy Clean Vinyl Surfaces
- Removable Bowl with Lid & Carry Handle
- Height Adjustable
- Raised Padded Armrests
- Fold Away, Padded Seat
- Non-Slip Tips
- Toilet Seat to Floor Height 435 - 585 mm

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Toileting

Category: Prescribed

Click [here](#) for more information

Warranty Information

- 1 Year



NDIS

Item Code:

05_091203048_0103_1_2

Toilet Aids – Composite:

05_414_0103_1_2

Click [here](#) for more information

Compliance

AS/NZS 3973:2009

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Seat Depth	Seat Width	Seat To Floor Height
BEC045115	Aspire Deluxe Bedside Commode -	600 mm	600 mm	865 - 1015	160 kg	12 kg	Champagne	450 mm	460 mm	410 - 565 mm

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Seat Depth	Seat Width	Seat To Floor Height
	Champagne			mm			Vinyl			
BEC045116	Aspire Deluxe Bedside Commode - Slate	600 mm	600 mm	865 - 1015 mm	160 kg	12 kg	Slate Vinyl	450 mm	460 mm	410 - 565 mm