

Cuddle Bed Configuration



- Cuddle Beds allow for hospital stays to feel like home, creating comfort for both patients and loved ones.
- [Benefits of co-sleeping in a hospital setting](#)
- Boasts a 400 kg SWL and electronic height adjustment range of 390 mm to 760 mm
- Electronic High/Low with backrest, knee bend and Trendelenburg
- Anti-jam backrest & knee bend motor
- 10-button soft touch pendant handset with auto profiling
- Two single bed join together and work independently or in unison
- OH&S friendly and aesthetically pleasing radius corners
- Simple connections between bed frame and electrics
- Backrest quick release CPR Manual Lever
- 4 x 100 mm central and directional lockable castors
- Easily transfer bed between rooms
- Control bed movement from either the left or right side of the bed
- Includes: Backup battery, IV pole & bolsters
- Built-in Foot Extension: 150 - 300 mm
- Hoist clearance: 150 mm
- Removable moisture-resistant MDF head & foot panels
- Mattress retainer instead of foot panel
- Perforated sheet metal mattress deck for easy cleaning
- Enhanced infection control & mattress breathability
- Epoxy coated & anti-rust treated steel frame
- Teflon-coated pivot bushes for smooth & quiet operation
- Safe low-voltage power supply

✓ Contract approved

Support @ Home

Function: Beds and bed equipment

Category: Prescribed

Click [here](#) for more information

Warranty Information

- Steel Construction - 15 years
- Electronics - 3 years
- Castors - 3 years

Code	Display Name	Total Length	Total Width	Safe Working Load
BEB045106	Cuddle Bed Configuration – 400kg	2060mm	1940mm	400kg