

## Jay - Zip (Paediatric)



- Dual layered contoured foam base
- Reticulated foam in the outer cover promotes air movement
- Soft and stretchable outer cover
- All models include: 2 way stretch cover

### Warranty Information

- 2 Years Cushion
- 6 Months Cover

**JAY®**



Contract approved

### Support @ Home

**Function:** Manage tissue integrity

**Category:** Prescribed

Click [here](#) for more information

### NDIS

**Item Code:**

05\_181006190\_0105\_1\_2

**Pressure Care Cushion –**

**Composite:**05\_429\_0103\_1\_2

Click [here](#) for more information

### Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material	Pressure Care Rating
PTC585184	JAY ZIP - 10 x 10in	10 in	10 in	2.75 in	75 kg	Foam	Medium Risk
PTC585193	JAY ZIP - 12 x 12in	12 in	12 in	2.75 in	75 kg	Foam	Medium Risk
PTC585185	JAY ZIP - 12 x 16in	16 in	12 in	2.75 in	75 kg	Foam	Medium Risk

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material	Pressure Care Rating
PTC585192	JAY ZIP - 14 x 14in	14 in	14 in	2.75 in	75 kg	Foam	Medium Risk
PTC585190	JAY ZIP - 14 x 16in	16 in	14 in	2.75 in	75 kg	Foam	Medium Risk
PTC585191	JAY ZIP - 15 x 15in	15 in	15 in	2.75 in	75 kg	Foam	Medium Risk