

## **Aspire Autoclavable Bed Pan**



Our Autoclavable Bed Pans are made from premium polypropylene. They will last for many years in any healthcare environment subject to the rigours of cleaning, sanitation and infection control requirements.

- Heavy duty and reusable
- Designed for use in healthcare environments such as surgeries, acute care, general ward, medical centres, and aged care facilities
- Autoclavable polypropylene that can be cleaned using autoclaves and sanitisers
- Pastel green colour for easy identification

## Warranty Information

• 1 Year



$\odot$	Contracted	DVA	product
			•

Ocontract approved

## ndis

**Item Code:** 

05\_098800044\_0103\_1\_2

**Toilet Aids – Composite:** 

05\_414\_0103\_1\_2

Click here for more information

## **Compliance**

AS/NZS 3973:2009

Code	Display Name	Total Length	Capacity
BTT151800	Bed Pan	360 mm	1000 ml