

## Liko Viking M Lifter



- Ideal for meeting the lifting needs in small spaces
- Informative control box displays number of lift cycles, battery level, lift overload and service indicator
- Optional quick-release hook system which means caregivers can easily and rapidly change a slingbar or install a scale
- Optional armrests help support patients during ambulation as well as making it easier for caregivers to manoeuvre the lift
- Three variable height settings accommodate a variety of different lifting situations
- Intuitive and ergonomically-designed hand control incorporates indicators for low battery level and when the Safe Working Load of the lift is exceeded
- An extensive assortment of accessories is available including slings, slingbars and scales
- Lifting Height 420 - 1790 mm

### Warranty Information

- Frames: 10 Years
- Electrical Parts: 2 Years



Contract approved

### Support @ Home

**Function:** Lifting persons

**Category:** Prescribed

Click [here](#) for more information

### NDIS

**Item Code:**

05\_123603161\_0105\_1\_2

**Hoist –**

**Composite:**05\_420\_0105\_1\_2

Click [here](#) for more information

### Compliance

AS/NZS ISO 10535:2011

| Code      | Display Name                          | Total Length | Total Width   | Total Height   | Safe Working Load | Material  |
|-----------|---------------------------------------|--------------|---------------|----------------|-------------------|-----------|
| LSS393190 | Liko Viking Lifter - MEDIUM           | 1230 mm      | 690 - 1110 mm | 1340 - 2020 mm | 205 kg            | Aluminium |
| LSS393200 | Liko Viking - M - Leg Protector Set - | -            | -             | -              | -                 | -         |