

## MedTime SAFE Tamper Resistant Automatic Pill Dispenser



**e-pill | MedTime SAFE | Tamper Resistant Automatic Pill Dispenser** Helps Prevent Double Dosing and Keep Medication Safe. It's a tamper resistant automatic pill dispenser with four independent locks and can dispense up to 28 times per day. This portable dispenser is a locked automatic medication dispenser for home or institutional use for pain, psych, or drug rehabilitation medicines.



- Robust device (with multiple locks and metal housing)
- Allows for “take-home” opioid treatment and easy dose verification by caregivers.
- Security: The metal case blocks access to the automatic pill dispenser and medication tray, creating a highly tamper-resistant environment.
- Easy to use – tip the pillbox over when the alarm sounds to dispense meds
- Tamper-resistant for pain meds, psych meds, and CIII medications
- Automatic Pill Dispenser housed in locked metal box
- Four locks (eight keys – two per lock included)
- Dispense up to 28 times per day with sound and visual notification
- Complete medication management system included
- No monthly fees, 30-day money-back guarantee, and 1-year warranty
- Capacity: Each compartment holds 18 aspirin-sized tablets. Max 28 doses per day.
- Keys: Two sets (8 total) included.
- Battery life: 6–12 months (4 × AA included) with low battery indicator.

**Please note:** There are no key overrides or duplicate keys available. Do not lose your keys.

Four (4) Locks are included with your e-pill MedTime SAFE:

- Pill dispenser lock (2 keys)
- Metal pill box lid lock (2 keys)

- Padlock (2 keys)
- Location lock (2 keys)
- Spare trays available for pharmacy refilling (optional extra)
- Tamper Evident Labels available to assist pharmacy compliance with QCPP requirements

| Code      | Display Name                                           |
|-----------|--------------------------------------------------------|
| DLG275020 | MedTime SAFE Tamper Resistant Automatic Pill Dispenser |
| DLG275030 | e-pill MEDTIME Locked Automatic Pill Dispenser         |
| DLG275260 | Spare Tray - e-pill Medtime Pill Dispenser             |