

M1 Corpus



permobil

- **Less fuss, more adaptable:** A quality solution for those users who require a more traditional rehab style seating system with limited to no power seat functions. Designed with adaptable seating and positioning needs in mind, the M1 offers all the flexibility without the compromise.
- M1 Core package available - [View Package Here](#)
- Locking flip back arms with easy to release mechanism.
- Multi-angle adjustable swing-away legrests.
- FlexLink™ suspension and precision engineered CasterGlide™ technology to help ensure smooth, confident operation in a variety of driving scenarios.
- Designed with 1" diameter rigid tubular backrest frame allows mounting of virtually any seating and positioning solution.
- Transport friendly features including standard crash tested transport brackets, optional EZ fold backrest and a lower overall chair weight of 145 kg with batteries, expands available travel options



Contract approved

Support @ Home

Function: Powered wheelchairs

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_122306125_0105_1_2

Wheelchair Powered –

Composite:05_419_0105_1_2

Click [here](#) for more information

Warranty Information

- 1 year

Code	Display Name	Total Width	Total Height	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height	Maximum Speed	Range
PWC652575	M1	24 in	18 to 29 in	136 kg	145 kg Inc. batteries	35 - 56 cm	14 - 22 in	40 - 46 cm	8 km/h	25 km