

SitnStand Lift Assist - Classic



Upgrade your home and existing chair with our Lift Assist Classic, designed to enhance mobility and independence. This versatile device offers reliable assistance for individuals who require support when transitioning between sitting and standing positions at a fraction of the cost of a lift chair.

- **Sturdy Construction:** Built with durable materials for long-lasting performance and stability.
- **Simple Installation:** Easily install without the need for complicated tools or professional assistance.
- [Watch the informational video here.](#)
- **Smooth Operation:** Experience smooth and effortless transitions from sitting to standing with the push of a button.
- **Comfortable Design:** The lift assist is designed with user comfort in mind.
- **Compact Size:** Its compact size allows for easy placement in various areas of your home, maximising convenience.
- **Weight Allowance:** Up to 266lbs/121kg.
- **Rechargeable Battery:** Last for 4-7 days.

Warranty Information

- 2 years

✓ Contract approved

Support @ Home

Function: Changing body position

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Safe Working Load
LSS403710	SitnStand Lift Assist - Classic	121 kg