



COVE PRESSURE RELIEF CHAIR

The Aspire Cove chair is designed to improve postural, pressure care, and comfort outcomes for users requiring medium to high level support.



Cove Chair
Brochure



Legrest
Release Lever



Ergonomic Carer
Push Handle



Optional
Portering Footrest
CHP198985



Optional
Tray Table
CHP198980

3 YEAR WARRANTY
Frame

1 YEAR WARRANTY
Fabric



- **Adjustable Pressure Seat Cushion:** Hybrid air-foam cell with simple valve adjustment
- **Adjustable Seat Depth:** Forward/back backrest positioning
- **Infinite Positioning:** 40° tilt and independent leg raise, both gas-assisted for easy operation
- **Infection Control:** Medical grade upholstery with biostatic medical surfaces
- **Contoured Head Support:** Adjustable and neck support pillow
- **Tool-free Adjustment:** Fast setup without tools
- **Tracking Armrests:** Follow tilt positioning for user confidence during repositioning
- **Ergonomic Handles:** Multiple grip positions for carers
- **Anatomical Legrest:** Extends outward following knee movement
- **High-Performance Castors:** TENTE Levina wheels reduce effort
- **Pressure Cushion Compatible:** Works with ROHO and similar systems

Code	Description	Overall Width	Overall Height	Overall Length	Seat Width	Seat Depth	Seat Height*	Backrest Height
CHP198978	Medium	780mm	1270mm	860 - 1640mm	550mm	500 - 580mm	540mm	740mm
CHP198979	Large	830mm	1300mm	860 - 1690mm	600mm	500 - 580mm	540mm	740mm

*Seat-to-floor Height