




# POSTUREPRO CHAIR RANGE


The PosturePro range offers highly adjustable seating solutions designed to deliver optimal pressure care and postural support for individuals who sit for extended periods.




**Customisable Backrest with Integrated Laterals:** Provides side support to encourage an upright posture and reduce leaning, helping maintain proper alignment



**Tool-free Seat Adjustment:** Seat depth and width can be quickly adjusted without tools for a better fit



**Adjustable Armrests:** Height and depth adjustment allows correct arm positioning and better support



**Adjustable Footrest:** Depth and angle adjustment supports posture and sit-to-stand transfers



**Standing Aid Compatible:** Works with Aspire Star Mover and most Aspire patient lifters for safe transfers



**Multiple Positioning Options:** Tilt and recline functions assist with repositioning, relieve pressure and help prevent pressure injuries

# ADDITIONAL FEATURES OF POSTURE**PRO** PLUS


**Enhanced Postural Support:** Adjustable head, neck, and shoulder sections provide tailored support for users with complex postures, including kyphosis and contractures

**Adaptable to Progressive Conditions:** Designed to adjust over time to meet changing support needs

**Negative Angle Leg Rest:** Supports users with contractures or specific leg positioning requirements


**Pocketed Waterfall Backrest:** Height and lateral angles are adjustable, with padding that can be added or removed for customised support

## CLINICAL CONSIDERATIONS



### Improved postural positioning to support

- Pressure care management
- Users with low tone or postural deviations; achieved through full body support and use of tilt-in-space, recline, and leg elevation to counter gravity forces
- Breathing
- Feeding, drinking, and swallowing
- Communication
- Sight line for users with kyphotic postures
- Meaningful occupational activities and social engagement
- Relaxation and fatigue management



### Pressure redistribution and pressure care support

- Tilt and recline functions
- Breathable, multi-stretch fabric on all support surfaces
- A combination of immersive memory foam and air bladders
- Customisable support surfaces (head, back, seat, arm, leg, and feet) to promote optimal posture and minimise shear and pressure forces



### Carer Support and reduced manual handling

- Ability to change the user's posture without returning them to bed
- Use of chair functions to assist with positioning for different activities
- Large wheels to enable easier transport of user to different environments
- Compatibility for use with active lifters and transporter aids
- Tilt-in-space and recline to support positioning and use of a full hoists and ceiling hoists
- Gas strut for easy use of tilt and recline functions

**2** YEAR WARRANTY  
Frame

**1** YEAR WARRANTY  
Components and Fabric



Code	Description	Overall Depth	Overall Width	Overall Height	Seat Depth	Seat Width	Seat Height*	Backrest Height	Backrest Tilt Angle Range	Seat Tilt Angle Range	Leg Elevation Range
CHP231600	PosturePro	950	650 - 850	1300 - 1360	450 - 540	330 - 530	530 - 590	770	90 - 130°	-5 - +32°	0 - 35°
CHP231610	PosturePro Plus	950	650 - 888	1111 - 1400	450 - 540	330 - 530	530 - 590	777 - 900	90 - 130°	-5 - +32°	-5 - 30°

\*Seat-to-foot plate Height