

aspire™ ComfiMotion®

Range Selection Guide

Activ Care Ultra Low Adjustable Bed

Clinical Positioning

Split-frame Design

Memory Function

Lockable Castors

Ultra-low Height

Customisable Safe Entry/Exit Height

Underbed Clearance

Battery Backup

Night Assist

Illuminated Handset

Massage Function



Hi Lo



Back
Raise



Head
Tilt



Leg
Raise



Zero
Gravity



Seated



Trendelenburg



Care Variable Height Adjustable Bed

Clinical Positioning

Lockable Castors

Night Assist

Illuminated Handset

Enclosed Frame Design

Safe Lowering Technology

Box-on-box Design



Hi Lo



Back
Raise



Leg
Raise



Seated



Trendelenburg



Relax Fixed Height Adjustable Bed

Memory Function

Night Assist

Massage Function

Press and Hold Safety Handset

Box-on-box Design

German Engineered Motors

Personalised Positioning

USB Charge



Back
Raise



Head
Tilt



Leg
Raise



Create your ComfiMotion®

We believe that creating a supportive and luxurious sleep environment can significantly improve health and wellbeing. Our ComfiMotion Collection is designed to blend seamlessly into the home environment, whilst improving comfort and independence.

To build your ComfiMotion configuration, simply select your bed base, mattress, headboard and accessories.



Beds

Care Adjustable Bed - (LS)	BEB046720
Care Adjustable Bed - (KS)	BEB046730
Relax Adjustable Bed - (LS)	BEB046900
Relax Adjustable Bed - (KS)	BEB046910
Relax Adjustable Bed - (LD)	BEB046920
Relax Adjustable Bed - (Q)	BEB046930
Activ Care Adjustable Bed - (LS)	BEB046745
Activ Care Adjustable Bed - (KS)	BEB046755
Activ Care Adjustable Bed - (LD)	BEB046765
Activ Care Adjustable Bed - (Q)	BEB046775
Fixed Companion Bed - (LS)	BEB046800

Mattresses

Pocket Spring - (LS)	BEM469800
Pocket Spring - (KS)	BEM469810
Pocket Spring - (LD)	BEM469820
Pocket Spring - (Q)	BEM469830
Memory Foam - Classic (LS)	BEM469910
Memory Foam - Classic (KS)	BEM469920
Memory Foam - Classic (LD)	BEM469922
Memory Foam - Classic (Q)	BEM469924
Memory Foam - Pressure Care (LS)	BEM469930
Memory Foam - Pressure Care (KS)	BEM469940
Memory Foam - Pressure Care (LD)	BEM469942
Memory Foam - Pressure Care (Q)	BEM469944

Pillows

Breeze - Memory Foam Cooling Gel	BEP050100
Plush - Memory Foam Fibrefill	BEP050110

Headboards Freestanding

Allure Headboard - (LS)	BEA830120
Allure Headboard - (KS)	BEA830130
Allure Headboard - (LD)	BEA830132
Allure Headboard - (Q)	BEA830134
Allure Headboard - (K)	BEA830136
Chesterfield Headboard - (LS)	BEA830100
Chesterfield Headboard - (KS)	BEA830110
Chesterfield Headboard - (LD)	BEA830112
Chesterfield Headboard - (Q)	BEA830114
Chesterfield Headboard - (K)	BEA830116

Headboards and Footboards Fixed*

Allure Headboard - (LS)	BEA830150
Allure Headboard - (KS)	BEA830152
Allure Headboard - (LD)	BEA830154
Allure Headboard - (Q)	BEA830156
Footboard - (LS)	BEA830160
Footboard - (KS)	BEA830162
Footboard - (LD)	BEA830164
Footboard - (Q)	BEA830166

Accessories

Crook Handle Bed Stick - Care/Relax Bed	BEA830200
Crook Handle Bed Stick - Activ Care Bed	BEA830211
Bed Rail - Care/Relax Bed	BEA830210
Bed Rail - Activ Care Bed	BEA830212
Extendable Bed Rail	BEA830205
Large Accessibility Castors - Care Bed	MWT451080
Battery Backup - Care/Relax Bed	BEA830220
Split King Upgrade Kit - Activ Care Bed	BEA830230
Simplified Handset - Activ Care Bed	BEA830214
Partner Handset - Activ Care Bed	BEA830213

*Fixed Allure headboards and footboards are compatible with ComfiMotion ActiveCare beds