

NUWEIGH

CHAIR SCALE

CHR-691 #1533

FEATURES

- Adjustable arm rests allow for easy entry/exit from the chair
- 300kg capacity
- 100gram accuracy
- Strong fibreglass chair
- BMI function
- Tubular design with no sharp edges
- Lift up swivel footrests for comfort and easier patient transfer
- Level indicator for consistent accuracy
- Large 26mm LCD display
- Built in rechargeable battery
- All cabling is protected
- Ergonomic 930mm handle height
- 4 Swivel castor wheels with locks
- Hygienic easy to clean seat



FUNCTIONS



Large 26mm
LCD display



Sturdy swivel
footrests



Adjustable
armrests



Comfortable
fibreglass chair

CAPACITY

300kg x 100g

PRODUCT CODE

LSS377720

ALWAYS FOLLOW THE DIRECTIONS FOR USE

NWS
NUWEIGH

Distributed by

aidacare
HEALTHCARE EQUIPMENT